

BBS Revolution™

The Next Generation Bladder Scanner

Fully automated bladder scanner with a unique, user friendly design.

- Automatically detects male/female anatomy
- Automatically locates bladder without aiming
- Quickly provides results you can trust

Improved Technology + Reduced Total Cost of Ownership +
Decreased Infection Risk
= BETTER SOLUTION



BBS Revolution Clinical Benefits

- ...Saves nurses time and frustration from setup through cleanup
- ...Reduces patient discomfort due to probe design and quick test result
- ...Simplifies disinfection requirements via the wireless operation
- ...Increases compliance with infection control protocols
- ...Improves accuracy and uptime with self-calibration mode
- ...Minimizes user training through its elegantly simple design

Waging War Against Healthcare-Associated Infections

- ✓ Lower Costs
- ✓ Better Outcomes
- ✓ More satisfied patients and providers

dBMED_x

BBS Revolution™

The Next Generation Bladder Scanner

Wireless
Scanner

8 transducers,
180 degree view angle

Internal Battery,
up to 8 hours continuous use

It takes just one quick scan to get an accurate, reliable volume measurement.

Follow these simple steps:

1. Apply ample gel on the abdomen
2. Place the probe firmly against the lower abdomen
3. Press the button to take a reading

IT IS THAT EASY!



There is no need to select modes or aim the probe, it is ONE AND DONE™ with the BBS Revolution™.

CONTACT US FOR A DEMO

- www.dBMEDx.com
- 720-515-6630
- info@dbmedx.com

dB MED X™